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AUTOMATIC ARM BLOOD PRESSURE MONITOR

USER MANUAL

Thank you for choosing our product. We hope you will enjoy using the appliance.

IMPORTANT INFORMATION

NORMAL BLOOD PRESSURE FLUCTUATION

All physical activity, excitement, stress, eating, drinking, smoking, body posture and many other activities or factors (including taking a blood pressure measurement) will influence blood pressure value. Because of this, it is mostly unusual to obtain identical multiple blood pressure readings.

Blood pressure fluctuates continually ----- day and night. The highest value usually appears in the daytime and the lowest one usually at midnight. Typically, the value begins to increase at around 3:00 AM, and reaches to highest level in the daytime while most people are awake and active.

Considering the above information, it is recommended that you measure your blood pressure at approximately the same time each day.

Too frequent measurements may cause injury due to blood flow interference, please always relax a minimum of 1 to 1.5 minutes between measurements to allow the blood circulation in your arm to recover. It is rare that you obtain identical blood pressure readings each time.

CONTENTS AND DISPLAY INDICATORS

- 1 Button MEM
- 2 Button START
- 3 Battery compartment
- 4 LCD display
- 5 Cuff
- 6 Low battery indicator
- 7 Systolic pressure
- 8 Diastolic pressure / pulse rate display (alternating)
- 9 Date / Time display (alternating)
- 10 Memory indicator
- 11 Irregular heartbeat symbol
- 12 Blood pressure level classification Indicator
- 13 "Ready to inflate" indicator

INTENDED USE

The fully automatic electronic sphygmomanometer is for use at home. It is a non-invasive blood pressure measurement system intended to measure the diastolic and systolic blood pressures and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm. The cuff circumference is limited to 22 cm - 30 cm.

CONTRAINDICATION

It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

PRODUCT DESCRIPTION

Based on Oscillometric methodology and silicon integrated pressure sensor, blood pressure and pulse rate can be measured automatically and non-invasively. The LCD display will show blood pressure and pulse rate. The most recent 2x60 measurements can be stored in the memory with date and time stamp. The monitor can also show the average reading of the last three measurements.

The Electronic Sphygmomanometers corresponds to the below standards: IEC 60601- 1:2005/EN 60601-1:2006/AC:2010 (Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance), IEC60601-1-2:2007/EN 60601-1-2:2007 /AC:2010 (Medical electrical equipment -- Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility – Requirements and tests), IEC 80601-2- 30:2009+Cor.2010/ EN 80601-2-30:2010 (Medical electrical equipment –Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) EN 1060-1: 1995 + A1: 2002 + A2: 2009 (Non-invasive sphygmomanometers - Part 1: General requirements), EN 1060-3: 1997 + A1: 2005 + A2: 2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems).

SPECIFICATIONS

1. Product name: Blood Pressure Monitor
2. Model: JC-1321 (KD-558)
3. Classification: Internally powered, Type BF applied part, IPX0, No AP or APG, Continuous operation
4. Device size: Ca. 138 mm x 48 mm x 98 mm
5. Cuff circumference: 22 cm-30 cm
6. Weight: Approx. 211 g (batteries and cuff excluded)
7. Measuring method: Oscillometric method, automatic inflation and measurement
8. Memory volume: 2 users x 60 measurements with time and date stamp
9. Power source: Batteries: 4 x 1.5 V  type AA
10. Measurement range:
 - Cuff pressure: 0 – 300 mm Hg
 - Systolic: 60 – 260 mm Hg
 - Diastolic: 40 – 199 mm Hg
 - Pulse rate: 40 – 180 beats/minute
11. Accuracy:
 - Pressure: 3 mm Hg
 - Pulse rate: 5 %
12. Environmental temperature for operation: +10 °C to +40 °C (50 °F to 104 °F)
13. Environmental humidity for operation: ≤ 90 % RH
14. Environmental temperature for storage and transport: -20 ° to 55 °C (-4 ° to 131 °F)
15. Environmental humidity for storage and transport: ≤ 90 % RH
16. Environmental pressure: 80 kPa to 105 kPa

17. Battery life: Approx 500 measurements
18. A list of all components belonging to the pressure measuring system, including accessories: Pump, valve, LCD, cuff, sensor

Note: These specifications are subject to change without notice.

NOTICES

1. Read all of the information in the operation guide and any other literature in the box before operating the unit.
 2. Stay still, calm and rest for 5 minutes before blood pressure measurement.
 3. The cuff should be placed at the same level as your heart.
 4. During measurement, neither speak nor move your body and arm.
 5. Measuring on same arm each time.
 6. Please always relax at least 1 or 1.5 minutes between measurements to allow the blood circulation in your arm to recover. Prolonged over-inflation (cuff pressure exceed 300 mm Hg or maintained above 15 mm Hg for longer than 3 minutes) of the bladder may cause ecchymoma of your arm.
 7. Consult your physician if you have any doubt about below cases:
 - 1) The application of the cuff over a wound or inflammation diseases;
 - 2) The application of the cuff on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present;
 - 3) The application of the cuff on the arm on the side of a mastectomy;
 - 4) Simultaneously used with other monitoring medical equipments on the same limb;
 - 5) Need to check the blood circulation of the user.
 8. ⚠ This Electronic Sphygmomanometer is designed for adults and should never be used for infants or young children. Consult your physician or other health care professionals before use on older children.
 9. Do not use this unit in a moving vehicle; this may result in erroneous measurement.
 10. Blood pressure measurements determined by this monitor are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard Institute, Electronic or automated sphygmomanometers.
 11. If an irregular heartbeat (IHB) from common arrhythmias is detected in the procedure of blood pressure measurement, this sign (♥) will be displayed. Under this condition, the Electronic Sphygmomanometers can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment.
- There are 2 conditions under which the signal of IHB will be displayed:
- 1) The coefficient of variation (CV) of the pulses is $> 25\%$.
 - 2) The deviation of the following pulse period is ≥ 0.14 s, and the number of such pulses amounts to more than 53 % of the total number of measured pulses.
12. Please do not use the cuff other than supplied by the manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.
 13. ⚠ The monitor might not meet its performance specifications or cause a safety hazard if stored or used outside the temperature and humidity ranges specified in the specifications.
 14. ⚠ Please do not share the cuff with other contagious persons to avoid cross-infection.
 15. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television

reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the device to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

16. This blood pressure monitor is verified by auscultatory method. It is recommended that you check annex B of ANSI/AAMI SP-10:2002+A1:2003+A2:2006 for details of verification method if you need.

NOTE: Do not forget: while self-measurement represents control, it does not constitute a proper diagnosis or treatment. Any unusual values should always be discussed with your doctor.

SETUP AND OPERATING PROCEDURES

1. BATTERY LOADING

- Open battery cover at the back of the monitor.
- Load four batteries of the type "AA". Please pay attention to polarity.
- Close the battery cover.

When the LCD shows the battery symbol , replace all batteries with new ones.

Rechargeable batteries are not suitable for this monitor.

Remove the batteries if the monitor will not be used for a month or more to avoid relevant damage of battery leakage.


 Do not let battery fluid get into your eyes. If it should get in your eyes, immediately rinse with plenty of clean water and contact a physician.

The monitor, the batteries and the cuff, must be disposed of according to local regulations at the end of their usage.

2. CLOCK AND DATE ADJUSTMENT

- Once you install the battery or turn off the monitor, it will enter Clock Mode, and LCD will display time and date by turns. See Fig. 2 and 2-1.

- While the monitor is in Clock Mode, pressing the "START" button and the "MEM" button simultaneously; a beep is heard and the month will flash first. See Fig. 2-2. Press the button "START" repeatedly; the day, the hour and the minute will flash in turn. While a number is flashing, press the button "MEM" to increase the number. Hold the button "MEM" depressed, the number will increase rapidly.

- You can turn off the monitor by pressing the "START" button when the minutes are flashing, then the time and the date are confirmed.


- The monitor will turn off automatically after 1 minute of no operation; with the time and date unchanged.

- After replacing the batteries, you should set the time and the date again.

3. CONNECTING THE CUFF TO THE MONITOR

Insert the cuff tubing connector into the socket in the left side of the monitor. Make certain that the connector is completely inserted to avoid air leakage during blood pressure measurements.

⚠ Avoid compression or restriction of the connection tubing during measurement, as this may cause an inflation error or harmful injury due to continuous cuff pressure.


4. APPLYING THE CUFF

a. If required, pull the cuff end through the medial loop (the cuff is packaged like this already).

b. Place the cuff around a bare arm 1 to 2 cm above the elbow joint.

c. While sitting, place your hand with the palm facing up onto a flat surface like a table or similar. Position the air tube in the middle of your arm in line with your middle finger.

d. Tighten the cuff and close the Velcro fastener. The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger between your arm and the cuff.

Note:

Please refer to the cuff circumference range in "SPECIFICATIONS" to make sure that the appropriate cuff is used.

1. Measuring on same arm each time.

2. Do not move your arm, body, or the monitor and do not move the rubber tube during measurement.

3. Stay quiet, calm for 5 minutes before blood pressure measurement.

4. Please keep the cuff clean. If the cuff becomes dirty, remove it from the monitor and clear it by hand in a mild detergent, then rinse it thoroughly in cold water. Never dry the cuff in clothes dryer or iron it. Clean the cuff after the usage of every 200 times is recommended.

5. BODY POSTURE DURING MEASUREMENT

Sitting Comfortably Measurement

a. Be seated with your feet flat on the floor, and don't cross your legs.

b. Place palm upside in front of you on a flat surface such as table.

c. The middle of the cuff should be at the level of the right atrium of the heart.


Lying Down Measurement

a. Lie on your back.

b. Place your arm straight along your side with your palm upside.

c. The cuff should be placed at the same level as your heart.

6. TAKING YOUR BLOOD PRESSURE READING

a. After applying the cuff and with your body in a comfortable position, press the "START" button. A beep is heard and all display characters are shown for self test. See Fig. 6. Please contact the service centre if a segment is missing.

b. Then the current user number is shown (U1 or U2) is flashing. See Fig. 6-1. Press the "MEM" button to change to the other user. See Fig. 6-2. Confirm your selection by pressing the "START" button. The current bank will also be confirmed automatically after 5 seconds with no operation.

6

6.1

6.2

- c. After selecting the memory user, the monitor starts to seek zero pressure. See Fig. 6-3.
- d. The monitor inflates the cuff until sufficient pressure has built up for a measurement. Then the monitor slowly releases air from the cuff and carries out the measurement. Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen separately. Irregular heartbeat symbol (if any) will blink. See Fig. 6-4 & 6-5. The result will be automatically stored in the current memory bank.

- e. After measurement, the monitor will turn off automatically after 1 minute of no operation. You can also press the "START" button to turn off the monitor manually.
- f. During measurement, you can press the "START" button to turn off the monitor manually.

Note: Please consult a health care professional for interpretation of pressure measurements.

7. DISPLAYING STORED RESULTS

This device is capable of memorizing data for two users: up to 60 readings for each user. U1 (user 1) - U2 (user 2) After each reading, the JOYCARE blood pressure monitor stores in memory all data from user 1 or 2 depending on the initial setting.


- a. After the measuring, you can review the measurements in the current memory bank by pressing the button "MEM". Now the LCD displays the amount of the results in the current bank. See Fig. 7.

- b. You can also press the "MEM" button in Clock Mode to display the stored results. The current memory bank will flash and the amount of results in this bank will be displayed. See Fig. 7-1. Press the "START" button to change to the other bank. See Fig. 7-2. Confirm your selection by pressing the "MEM" button. The current bank will also be confirmed automatically after 5 seconds with no operation.
- c. After selecting the memory bank, the LCD will display the average value of the last three results in this bank. See Fig. 7-3 & 7-4. If no results are stored, the LCD will show dashes as shown in Fig. 7-5.

- d. When the average is displayed and you press the "MEM" button, the most recent result will be displayed. See Fig. 7-6. Then the blood pressure and the pulse rate will be shown individually. Possibly the irregular heartbeat symbol may flash. See Fig. 7-7 & 7-8. Press the "MEM" button again to display the next result. See Fig. 7-9. In this way, repeatedly pressing the "MEM" button displays the respective results measured previously.

e. When displaying the stored results, the monitor will turn off automatically after 1 minute of no operation. You can also press the “START” button to turn off the monitor manually.

8. DELETING MEASUREMENTS FROM THE MEMORY


When any result (except average reading of the last three results) is displayed and you keep pressing the “MEM” button for three seconds, all results in the current memory bank will be deleted after three “beeps”. The LCD will show Fig. 8; pressing the “MEM” or the “START” button will turn off the monitor.

9. ASSESSING HIGH BLOOD PRESSURE FOR ADULTS

The following guidelines for assessing high blood pressure (without regard to age or gender) have been established by the World Health Organization (WHO). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) also must be taken into consideration. Consult with your physician for accurate assessment, and never change your treatment by yourself.

Blood pressure classification for adults

Blood pressure classification	SYS (mm HG)	DIA (mm HG)
Optimal	< 120	< 80
Normal	120-129	80-84
high - normal	130-139	85-89
Hypertension level 1	140-159	90-99
Hypertension level 2	160-176	100-109
Hypertension level 23	≥ 180	≥ 110

Definition and classification of the blood pressure values according to WHO/ISH

10. TROUBLESHOOTING (1)

Problem	Possible cause	Solution
LCD Display shows abnormal result	The cuff position was not correct or it was not properly tightened	Apply the cuff correctly and try again

Problem	Possible cause	Solution
LCD Display shows abnormal result	Body posture was not correct during testing	Review the section "BODY POSTURE DURING MEASUREMENT" section of the instruction and try again.
	Speaking, arm or body movement, angry, excited or nervous during testing	Re-test when calm and without speaking or moving during the test
	Irregular heartbeat (arrhythmia)	It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

11. TROUBLESHOOTING (2)

Problem	Possible cause	Solution
LCD shows low battery symbol 	Low battery	Low battery
LCD shows "Er 0"	Pressure system is unstable before measurement	Don't move and try again.
LCD shows "Er 1"	Fail to detect systolic pressure	
LCD shows "Er 2"	Fail to detect diastolic pressure	
LCD shows "Er 3"	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again
LCD shows "Er 4"	Pneumatic system leakage or cuff is too loose during inflation	
LCD shows "Er 5"	Cuff pressure above 300 mm Hg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor
LCD shows "Er 6"	More than 3 minutes with cuff pressure above 15 mm Hg	
LCD shows "Er 7"	EEPROM accessing error	
LCD shows "Er 8"	Device parameter checking error	
LCD shows "Er A"	Pressure sensor parameter error	
No response when you press a button or load batteries.	Incorrect operation or strong electromagnetic interference.	Take out batteries for five minutes, and then reinstall all batteries.

MAINTENANCE

1.  Do not drop this monitor nor subject it to strong impacts.
2.  Avoid high temperatures and direct sunlight. Do not immerse the monitor in water as this will result in damage to the monitor.
3. If this monitor was stored at a temperature near the freezing point, allow it to come to room temperature before use.
4.  Do not attempt to disassemble this monitor.
5. If you do not use the monitor for a long time, please remove the batteries.
6. It is recommended the performance should be checked every 2 years or after repair. Please contact the service centre.
7. Clean the monitor with a dry, soft cloth or a soft cloth squeezed well after moistened with water, diluted disinfectant alcohol, or diluted soap suds.
8. No component in the monitor can be serviced by the user. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairably can be supplied.
9. The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years, and the cuff can maintain the performance characteristics for a minimum of 1000 measurements.
10. It is recommended that the cuff should be disinfected 2 times every week if needed (for example, in a hospital or in a clinic). Wipe the inner side (the side that contacts skin) of the cuff with a soft cloth moistened with ethyl alcohol (75-90%) and squeezed, then dry the cuff by airing.

EXPLANATION OF SYMBOLS ON UNIT

Symbol for "THE OPERATION GUIDE MUST BE READ". (The sign background colour: blue. The sign graphical symbol: white)

Symbol for "WARNING"

Symbol for "TYPE BF APPLIED PARTS" (the cuff is a type BF applied part)

Symbol for "ENVIRONMENT PROTECTION" - Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist.

Check with your local Authority or retailer for recycling advice.

Symbol for "MANUFACTURER"

Symbol for "DATE OF MANUFACTURE"

Symbol for "EUROPEAN REPRESENTATION"

SN Symbol for "SERIAL NUMBER"

Symbol for "KEEP DRY"

CE 0197

JC-1321 (KD-558) has been manufactured in compliance with the European directives 93/42/EEC 0197

MADE IN CHINA Imported and distributed by JOYCARE brand, by JOYCARE S.p.A. socio unico, Via Fabio Massimo, 45 Roma- Italy

Manufacture: ANDON HEALTH CO., LTD. Add: No. 3 JinPing Street, Ya An Road, Nankai District, Tianjin 300190, China.

Lotus Global Co., Ltd. 1 Four Seasons Terrace West Drayton, Middlesex, London, UB7 9GG, United Kingdom Tel: +0044-20-75868010 Fax: +0044-20-79006187

DISPOSAL The device (including its removable parts and accessories) must not be disposed of together with municipal waste at the end of its life, but in compliance with European Directive 2002/96/EC. Since it must be handled separately from household waste, either carry it to a separately-collected waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new device with the same purpose. Any infringement will be severely prosecuted. The batteries used in this device must be disposed of in the special bins at the end of their life.

All the instructions and drawings are based on the latest information available at the time of going to print and may be subject to change.

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ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration - electromagnetic emissions		
The JC-1321 [KD-558] is intended for use in the electromagnetic environment specified below. The customer or the user of the JC-1321 [KD-558] should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The JC-1321 [KD-558] uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The JC-1321 [KD-558] is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2
For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic immunity			
The JC-1321 [KD-558] is intended for use in the electromagnetic environment specified below. The customer or the user of the JC-1321 [KD-558] should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity			
The JC-1321 [KD-558] is intended for use in the electromagnetic environment specified below. The customer or the user of the JC-1321 [KD-558] should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the JC-1321 [KD-558], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ 800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the JC-1321 [KD-558] is used exceeds the applicable RF compliance level above, the JC-1321 [KD-558] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the JC-1321 [KD-558].

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 4
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the JC-1321 [KD-558]			
The JC-1321 [KD-558] is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the JC-1321 [KD-558] can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the JC-1321 [KD-558] as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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GB: WARRANTY CARD

Joycare S.p.A., a company with a sole shareholder, offers a period of 2 year warranty after the date of purchase. The warranty does not cover damage caused by accidental breakage, transportation, incorrect maintenance or cleaning, misuse, tampering or repairs carried out by non authorised third parties. The following is also excluded: damage caused by incorrect product installation, and damage caused by wear and tear, obviously including parts subject to wear and tear. In order to obtain the replacement/repair of products under warranty in cases other than those mentioned above, but which are instead caused by manufacturing defects, please contact your local importer/distributor. Goods that are not held to be in good condition because of the aforesaid causes must be accompanied by this certificate, duly filled-in, as well as a receipt for the goods (invoices and transport documents) and purchasing documentation.

Please keep the original packaging of the product.

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