

JC-117

JOYCARE®

CE 0434  
(dir. 93/42/EEC)

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GB

### INSTRUCTIONS FOR USE AND MAINTENANCE

JC-117 is a unit for the treatment and relieve of respiratory ailments. It is designed principally for domestic use.

JC-117 is equipped with a powerful piston compressor that does not need lubrication. The unit includes an air filter to be periodically replaced (every 70 applications). It is built in compliance with the current European Manufacturing and Safety Standards for electromedical equipments (EN 60601-1; EN13544-1). JC-117 is supplied with the following accessories: nebulizing bulb with mouthpiece, nasal prong, (paediatric mask), adult mask, connection tubing, filter (3).

### SAFEGUARDS

- Never use an extension cord or electrical adapter. Keep the cord away from heated surfaces.

- Do not handle the plug with wet hands, and not use the product while bathing. Never submerge the unit in to water; should this occur, unplug immediately. Do not reach for a product that has fallen into water; unplug it first. After removing the product from the water, do not use it again. (Return it immediately to an authorized service station). Always unplug the power cord immediately after use. Service should be performed only by authorized repair voids the warranty. Follow the regulations in force for disposal. In view of their length, the power cord and connection hose could constitute a strangulation hazard.

## INSTRUCTIONS FOR USE AND CLEANING

Plug the power cord into the electrical outlet corresponding to the voltage of the unit.

To use the ampul and its accessories, consult the additional sheet.

When the application has finished, turn the unit off and carefully clean the accessories.

Do not wash the unit under running water or do not submerge it in water; wipe it only with a damp cloth using mild detergents. Do not use abrasive cleaners.

Replacing the spare filter: lift the filter support lid (3) by levering with a flat-ended screwdriver as illustrated in figure 1a; detach the filter to be replaced and insert the spare filter as illustrated in figure 1b.

Finally, insert the filter support lid (3) into its seat, making it adhere well to the surface.

In compliance with EMC standard (EN60601-1-2) Unit case not water – proof.

The unit should not be used in the presence of anesthetic mixture inflammable with air, oxygen or nitrous oxide.

## IMPORTANT SAFEGUARDS

This is a medical device and is to be used only upon medical prescription and to be operated in accordance with this manual. It is important that the patient reads and understands the information provided for use of the unit. Please contact your retailer for any question.

Children and invalids should always use the unit under strict supervision of an adult who had read this manual.

## MICROBIC CONTAMINATION

In the presence of pathologies with the risk of infection and microbic contamination, a personal use of the accessories and the nebulizing bulb is suggested (always consult your physician).

### DISPOSAL

DISPOSAL PROCEDURE (Dir. 2012/19/Ue-WEEE) The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

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## GB GENERAL INFORMATION

Designed and created in accordance with the latest standards, NEBJET is an Medical Device (Risk Class IIa - MDD 93/42/EEC) device highly efficient in medical therapies administered by means of aerosol nebulization.

### Safety Warnings

- A single-patient use is recommended
- The device should be used only with compressor nebulizers for aerosoltherapy that comply with standards in force and whose use is indicated in the instruction manual
- Device not suitable for anaesthesia and lung ventilation
- The device should be used only with original accessories shown in this manual
- Information concerning connection and use with devices for aerosoltherapy are shown in the instruction manual of the devices. Always observe instructions and safety warnings as indicated in the instruction manuals of the devices for aerosoltherapy
- Use the device and its accessories as your doctor's recommendations
- Use medications prescribed by your doctor according to his indications
- Do not use the device for other use than Nebulizer for Aerosoltherapy. Manufacturer will not be held responsible for inappropriate use of the device
- The device does not come in sterile pack; always carry out cleaning and sterilization operations before and after each treatment. When using the device follow ordinary hygiene and cleaning precautions, (wash your hands and carry out cleaning operations as described in the specific section of this manual)
- The device has small components which might be removed and easily swallowed. Use by minors and disabled people require presence of an adult with his facilities. Do not leave the device unattended in places easily accessible by minors and disabled people
- Always store the device in a dry and clean place far from light, heat and weather
- Dispose of device as provided by standards in force
- In presence of pathologies with microbial contamination and infection hazard a singlepatient use of accessories and nebulizer is recommended (always consult your doctor)

How to prepare and to operate the device

The device is not sterile. Before use carry out cleaning and

disinfection operations as described in the specific section to pour the medication rotate the nebulizer top counter clockwise and remove it. For the quantity of medication prescribed by your doctor into the nebulizer tank. **IMPORTANT:** In case you do so, empty the nebulizer, clean it thoroughly and proceed with new filling. When medication has been poured, connect the nebulizer top onto the tank, rotate clockwise and insert the accessory prescribed by your doctor into the nebulizer top. Firmly connect one end of the air tube to the nebulizer and the other end to the air outlet on the compressor nebulizer for aerosoltherapy. Start the compressor for treatment to begin.

**IMPORTANT: NEVER INHALE IN HORIZONTAL POSITION. NEVER BEND THE NEBULIZER OVER 60°.**

### Cleaning instruction

After each treatment clean thoroughly each component of the nebulizer removing medication residuals and possible impurities. Clean each component as described here below except the air tube which -in case of treatment to a new patient or in case of impurities should be replaced. Wash all components (except air tube) under tap warm water (max 60 °C ) for about five minutes adding if necessary a small quantity of detergent following dosage and use limitations as provided by detergent manufacturer. Rinse thoroughly making sure that all deposits are washed away and let dry. All nebulizer components can be disinfected (except air tube) with chemical disinfectants following dosage and use limitations as provided by disinfectant manufacturer . They are usually available at pharmacy shops (Milton®, Amuchina® etc.). All nebulizer components can be heat steam sterilized (except air tube and masks) to 121°C (20 min.) or 134°C (7 min.) (EN554/ISO11134). The Sterilisation packaging must be conform to EN868/ISO11607 and be suitable for steam sterilisation. After sterilization always let all nebulizer components cool down to ambient temperature before use. Do not repeat sterilization cycle when nebulizer components are still warm.

### Nebulizer life

It is advisable replace NEBJET nebulizer after some 100 to 120 treatments on single patient or after about 20 sterilization cycles

Storage and transportation conditions		
Storage and transportation temperature:	- 20°C (-4°F) ÷ +70°C (158°F)	
Storage and transportation humidity percentage:	5% ÷ 95% RH	
WORKING TECHNICAL FEATURES		
Operating Minimum flow	5 l/min	+/- 10%
Operating pressure	60 kPa (0.6 bar, 8.7psi)	+/- 10%
Operating Max flow	7.5 l/min	+/- 10%
Operating pressure	150 kPa (1.5 bar, 21.75psi)	+/- 10%
Min. capacity	2 ml	
Max. capacity	2 ml	

Read carefully instructions for use

NEBJET conforms with the European directive 93/42/CEE dealing 0434 with medical devices, and to all the applicable European regulations.

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### GB: WARRANTY CARD

Joycare S.p.A., a company with a sole shareholder, offers a period of 2 year warranty after the date of purchase. The warranty does not cover damage caused by accidental breakage, transportation, incorrect maintenance or cleaning, misuse, tampering or repairs carried out by non authorised third parties. The following is also excluded: damage caused by incorrect product installation, and damage caused by wear and tear, obviously including parts subject to wear and tear. In order to obtain the replacement/repair of products under warranty in cases other than those mentioned above, but which are instead caused by manufacturing defects, please contact your local importer/distributor. Goods that are not held to be in good condition because of the aforesaid causes must be accompanied by this certificate, duly filled-in, as well as a receipt for the goods (invoices and transport documents) and purchasing documentation.

**Please keep the original packaging of the product.**

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