modello / model / modelo  
JC-230

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## EN User Manual

Read carefully the instruction before using the thermometer

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**GB: WARRANTY CARD**  
Joycare S.p.A., a company with a sole shareholder, offers a period of 2 year warranty after the date of purchase. The warranty does not cover damage caused by accidental breakage, transportation, incorrect maintenance or cleaning, misuse, tampering or repairs carried out by non authorised third parties. The following is also excluded: damage caused by incorrect product installation, and damage caused by wear and tear, obviously including parts subject to wear and tear. In order to obtain the replacement/repair of products under warranty in cases other than those mentioned above, but which are instead caused by manufacturing defects, please contact your local importer/distributor. Goods that are not held to be in good condition because of the abovementioned causes must be accompanied by this certificate, duly filled-in, as well as a receipt for the goods (invoices and transport documents) and purchasing documentation. Please keep the original packaging of the product.

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9. **1. INTRODUCTION**

Dear Client, thank you for buying FIDATI, the evolution of Thermofocus® - the world's first non-contact thermometer - and of VisiFocus®. Fidati is truly easy to use. It is capable of measuring a child or an adult's temperature without ever coming into contact with the skin: just bring it close to the forehead, at the distance that the thermometer itself will tell.

The FIDATI thermometer does not need to be placed in any parts of your child's body. If your baby is sleeping, you can use FIDATI without waking him up; and, if the child is awake, FIDATI will not bother him.

10. **2. WARNINGS**

Read carefully before using FIDATI

11. **2.1 Precautions**

1. Use FIDATI in a draft-free room, at a steady temperature between 16 and 40°C (60.8 and 104°F).

2. If the device comes from a room having an ambient temperature different from that of the room in which you are taking the reading (from another room, from a drawer, etc.), it needs to be stabilized. Do the MQCS before using it (see par. #6). As alternative, wait at least 5 minutes for the temperature to stabilize through the internal AQCS, (see par. #6) without touching the device.

3. Do not take a temperature reading if the person is sitting in a draft or if the subject has:

- been walking, running or exercising;

- come from another room that was ventilated or at a different temperature than the room where the thermometer is used;

- been wearing a cap, hat or scarf;

- been exposed to agents that could alter forehead temperature, e.g. shower, shampoo, hair-drier, sponging, etc.;

3. In all the above cases, wait a few minutes for the forehead temperature to stabilize.

4. Changing the reading point will bring to different results. Therefore, remember, always aim the projection on the same spot, precisely **at the center of the forehead** (midway between the top of the nose and the hairline) and keep the thermometer perpendicular to the forehead. **Do not take measurements on areas other than the center of the forehead**, except for the case at #4.2.

5. The temperature reading is taken in the area where the temperature is projected. It is of major importance to make sure that this area does not include eyebrows, hair or clothing. If necessary, brush away any hair from the forehead but remember, this must be done a few minutes beforehand or the temperature reading will be higher than the actual body temperature.

6. Avoid knocking and dropping it, and do not use it if damaged or if not functioning properly.

12. **2.2 Attention**

1. **Failure to observe the above-mentioned precautions (par. #2.1) may lead to very low or very high temperature readings, which cannot be attributed to product's malfunctions.**

2. The aiming lights meet the photo-biological safety requirements outlined in standard EN 62471. **No harm can be caused should the aiming lights be accidentally pointed in the eyes: the beams are harmless!**

3. The unit is a delicate measurement instrument and must not be used by little children. It is not a toy. Keep it out of the reach of children or persons with limited sensorimotor skills. Small parts can be ingested or inhaled.

4. Using this thermometer should not substitute medical consultation. Tell your physician what type of thermometer you are using and in which part of the body the temperature reading was taken.

5. If the child's fidgeting makes it difficult to take a correct reading, take first confidence with the device in order to approach the baby and take the temperature very quickly.

6. When taking your own temperature, use a mirror, or the black screen of a smartphone, or, if you have a child over six years of age, you can teach him how to use FIDATI.

7. Since it never comes into contact with the body, FIDATI does not require any "disposable" protection covers.

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05 - Italiano -

#### 4.3 Other readings

FIDATI can also be used to read the temperature of objects, food, liquids and other surfaces in the 1-55°C (33.8-131°F) temperature range. For example:  
 1. temperature of a baby's feeding bottle (fig. 7), food, bath, etc. In this case, mix the liquid well before taking the reading. With hot liquids or foods, take the reading quickly to prevent condensation from forming on the lens and wait 30 min. before taking another reading;  
 2. room temperature, by pointing the thermometer against a wall or a piece of furniture (see also par. #5);  
 3. temperature of a wound, inflammation, scar (for medical use).

**Proceed as you would do for a forehead reading (par. #4.1), but press the "home" button ( ); the projection will appear on the surface and the light of the "home" button will turn itself on.**  
 After the measurement, the projection of the temperature value will remain visible for about 12 seconds. During this time, should you find it difficult to see the temperature value, (as in the case of transparent, too dark or irregular surfaces, such as bathwater, soup) then, once released the button, you may address the projection on another light and opaque surface. You will need to respect the correct distance from the surface: indeed in this case, the thermometer is not taking a temperature, but only projecting the temperature value that has already been detected.

#### 5. MEMORY FUNCTION and AMBIENT TEMPERATURE

The "Mem" button lets you see the ambient temperature and call up the last 9 temperature readings (memory function).

- Aim the device at any uniform and clear surface.

- Press once the "Mem" button: the projection will show the letter "A" (Ambient) alternating with the ambient temperature value, and the light of the "Mem" button will turn itself on.

- Pressing again the "Mem" button, the projection will show the number 1 in alternation to the last temperature taken. The light of the button which had been used for that measurement ("faccia" or "caseta") will turn itself on. If you push the button again before the memory's light turns itself off, the number 2, alternating with the penultimate temperature taken will be projected, with the lightening of the used button will turn itself on. And so on.

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06 - Italiano -

#### 6. ROOM TEMPERATURE CALIBRATION

If, pressing the "face" ( ) button, the projection shows "CAL" [CAL].

Therefore, the thermometer will not let you take a measurement, which would not be correct. On the contrary it will activate an Automatic Quick Calibration System (AQCS). Now you have two possibilities:

1. Wait about 5 minutes without touching the thermometer, until the AQCS ends.

2. Perform a Manual Quick Calibration System (MQCS) as follows:

- press the "face" ( ) and "home" ( ) button simultaneously (fig. 8) as suggested by their flashing lights: the light of the "home" button will flash;

- within 10 seconds, focus the thermometer on an internal wall or wardrobe with uniform temperature and at a point between approximately between 80 and 150 cm (30 and 60 inches) from the floor and press the ( ) "home" button (fig. 9).

- once the right distance is reached (temperature values placed between the two dots fig. 4), release the button: the projection will flash slowly showing the room temperature.

To ensure a reliable temperature reading, do not focus the thermometer on an outside wall, window, source of heating or cooling (radiator, air conditioner, lamp, computer, surface in contact with the human body, etc.).

The thermometer is now ready to take a reading and will keep the MQCS for 15 minutes. The MQCS is possible only if the room temperature is in the 16-40°C (60.8-104°F) temperature range.

**Manual quick calibration (MQCS) can also be performed without the thermometer's suggestion if, for example, you need to move through rooms having different temperatures.**

#### 7. REPLACING the BATTERIES

• Set your thumb in the oval hollow on the back of the unit, press down and slide the battery hatch out as shown in figure 10.

• Remove the battery hatch.

• Remove the old batteries and dispose of them as required in the containers provided for this purpose.

• Insert 2 new AAA - LR03 batteries, preferably alkaline, carefully complying with the position indicated in their housing.

• To close the hatch, slide it in the opposite direction from which it was opened. After changing the batteries, let the thermometer stabilize for 20 minutes before taking a temperature reading, or run a MQCS (par. #6).

• Remove the batteries if you do not expect to use the thermometer for a long time.

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