## MANUALE D'ISTRUZIONE

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**INSTRUCTION MANUAL**

*Thanks a lot for having chosen a Joycare bathroom scale. This is a useful tool in order to check your weight, whenever you want.*

**TECHNICAL SPECIFICATIONS**

Max capacity: **180Kg/lb/st**

Graduation: 100 g / 0,2 lb/0.2lb (1lb $\geq$  20st)

GLASS platform

4 Strain Gauge accuracy sensors

Large and easy to read REVERSE display

Overloading indicator

Low battery indicator

Kg/lb selection

Power supply: 4xAAA batteries included

**OPERATING INSTRUCTION**

Select the weight unit kg/lb/st by the key placed near the battery compartment.

Put the scale on a hard, flat surface (avoiding carpets). The following scale turns on in response to pressure. Step on the scale; the display show your weight. The scale will switch off automatically after few seconds.

**To repeat the measurement**

Step down the scale and wait for display turning off; step on the scale again and measure your weight.

**Attention:** Stand still on the scale, distributing your weight evenly. Do not move.

**Battery replacement**

Open the battery case and insert 4 batteries paying attention to the polarity.

**MAINTENANCE AND CARE**

Please clean the unit with a damp and soft cloth avoiding any water seepage inside the scale. Do not use chemical detergents. Never deep the product into water. Never step on the scale platform with wet foot. Do not store the scale in an upright position when not in use; this may cause battery exhaustion. Avoid shaking, dropping or shocking the scale. This is a precision instrument and must be handled with care. Never try to repair the scale by yourself but always address to specialised centres. On the contrary, the warranty will be not valid anymore.

When the unit is not used for a long period of time, remove the battery, as matter of fact its leaking can damage the appliance. As quality of batteries varies from brand to brand, the product warranty does not include damages to the product caused by battery fluid leakage.

When the scale is off, never place or lay any objects on its platform. Always place and store the product away from heating sources, sunlight, dampness and sharp objects.

**In case of problems**

Check battery have been correctly fitted.

Check that you have selected the weigh unit.

Check whether the scale is on a flat, level floor and not touching against a wall.

Ensure whether you do not step on the scale until the display is turned off.

**Low battery indicator**

Whether the display shows "Lo" or it does not light up, batteries are flat. Please replace batteries.

**Overload indicator**

If "Err" is displayed, the scale has been overloaded (over 180 kg)

**WARNINGS**

Use this product only for its intended purpose, as described in this instruction manual such as electronic bathroom scale for domestic use only. Any other use should be considered improper and dangerous. Packaging components (plastic bags, cardboard, polystyrene etc.) must be kept out of the reach of children as they have the potential to be dangerous and their disposal should be carried out according to the current regulations. These scales should not be used to weigh objects or substances in commercial transactions, to produce medicines, to calculate tolls, tariffs, taxes, premiums, fines, remunerations, indemnities or fees of a similar kind determined by weight.

These scales are to be used for keeping one's own personal weight under control, and are not to be used for diagnostic purposes or for medical treatment.

Values which may seem unusual or outside the norm should always be discussed with your doctor.

**In order to avoid danger to the user of this device, the doses of any drug or therapy prescribed by one's own doctor should in no circumstances be modified on the basis of data calculated by this device.**

**NB: to avoid possible injury, do not step on the edge of the platform**

The product is made conforming with all the applicable European regulations.

**DISPOSAL**

The device (including its removable parts and accessories) must not be disposed of together with municipal waste at the end of its life, but in compliance with European Directive 2002/96/EC. Since it must be handled separately from household waste, either carry it to a separately-collected waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new device with

## **GB** JC-327 PERSONAL SCALE

the same purpose. Any infringement will be severely prosecuted. The batteries used in this device must be disposed of in the special bins at the end of their life.

**NOTE:** Specifications and designs are based on the latest information available at the time of printing and subject to change without notice.

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**GB: WARRANTY CARD**

Joycare S.p.A., a company with a sole shareholder, offers a period of 3 year warranty after the date of purchase. The warranty does not cover damage caused by accidental breakage, transportation, incorrect maintenance or cleaning, misuse, tampering or repairs carried out by non authorised third parties. The following is also excluded: damage caused by incorrect product installation, and damage caused by wear and tear, obviously including parts subject to wear and tear. In order to obtain the replacement/repair of products under warranty in cases other than those mentioned above, but which are instead caused by manufacturing defects, please contact your local importer/distributor. Goods that are not held to be in good condition because of the aforesaid causes must be accompanied by this certificate, duly filled-in, as well as a receipt for the goods (invoices and transport documents) and purchasing documentation.

**Please keep the original packaging of the product.**

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**JOYCARE S.p.A. socio unico**