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Thanks for having purchased a Joycare product.

Please for a correct usage of the product read carefully the following operating instructions.

It is advisable to preserve this manual for a future use.

USE AND REFERENCE STANDARD

This digital electric thermometer model is suitable to take body temperature for clinical and hospital use as well as for home use, as long as medical or paramedical staff are not required to be present at all times. The device was designed for oral, axillary or rectal use strictly in compliance with the instructions specified in this manual. Any use different from the use specified will be deemed improper use. The product complies with the following standards:- EN 60601-1 - Electromedical Appliances. Part 1: General Safety Regulations; - EN 60601-1-2 - Electromedical Appliances. Part 1: General Safety Regulations. Part 2 - Collateral Regulation: Electromagnetic Compatibility - Provisions and Tests; - UNI EN 12470-3 Clinical Thermometers - Compact Electric Thermometers Equipped with Maximum Temperature Device.

The digital electric thermometer model is a medical device belonging to Class IIa. Classification according to Clause 3.2, Section 3, Rule 10 of Directive 93/42/EEC.

SAFETY INSTRUCTIONS

Never bite, bend, drop or disassemble the thermometer. Keep the product out of children reach.

Avoid any contact with solvents. Do not expose the thermometer to direct sunlight, high temperature, moisture or

dirt. This thermometer has been designed exclusively for taking body temperature; any other use is improper and therefore dangerous.

Do not attempt to disassemble the unit, except for replacing flat battery

OPERATING INSTRUCTIONS

Before starting using the product, disinfect the probe with alcohol and then place the thermometer for taking the temperature. (before using the product, it is advisable to avoid your baby drinking hot or cold liquids and taking a bath or a shower).

Place the thermometer probe in the measuring place. Press the button and the display will show: 188.8°C for about 3 seconds, the last measurement

for about 3 seconds (if there is no previous measurements will be displayed \square L \square °C), then 37°C at the end the symbol °C (flashing for some seconds). Now the thermometer it's ready to measure your temperature. Once the measurement has ended, you hear a beep sound. Switch off the unit pressing the button.

Note : if the room temperature is above 32.0°C, the room temperature will display instead of the Lo

Oral usage: before taking the temperature, keep your mouth closed for 5 minutes. Place the probe under your tongue, close your mouth and switch the unit on. Keep your mouth well closed during the temperature reading, till the peak temperature has stabilised and you hear the beep alarm.

Note: normal temperature up to 37°C .

Axillary usage: wipe carefully your armpit with a dry towel, place the probe under the armpit and switch the unit on. Keep it in use up to 2 minutes, even if you hear the beep alarm.
Note: normal temperature up to 36.6°C.

Rectal usage: lubricate the probe with a water-soluble jelly. Do not use jelly like vaseline. Insert the probe into the rectum about 1 cm and switch the unit on. Wait till the temperature stabilises and you hear the beep alarm. Normal temperature: up to 37.6°C .

MEMORY FUNCTION

Thermometer stores the last temperature reading.

BATTERY REPLACEMENT

When the LCD displays the symbol “▼”, it would mean that battery needs to be replaced.

Use a screwdriver, in order to open the battery cover , and remove the old battery. Place a new LR41 or similar battery, respecting the polarity and close. When the unit is not used for a long period of time, please remove batteries, as matter of fact their liquid can damage the appliance.

WARNING: batteries must always be out of children reach and away from heating sources.

PRODUCT CLEANING

Clean always the thermometer, before or after each usage, with a soft cloth dumped with water.

Store the product into its plastic case. Never deep the unit into water or any other liquids.

TECHNICAL FEATURES

Measuring range	32.0 °C~42.9°C (90.0°F-109.9°F)
Accuracy	37.0°C-39.0°C ± 0.1°C; the rest ± 0.2°C;
Measuring time	15/20 seconds for oral and rectal use, about 30/40 for axillary use
Display	LCD
Memory	1
Power supply	Internally powered equipment LR41 or similar batter
Battery life	approx 1000 measurings
Usage conditions	Temperature: 10°C ~ +35°C Humidity: 15% ~ 80% RH
Storing, carriage and conservationconditions	Temperature: -10 ~ +60°C - Humidity25~90%RH

INFORMATION ON ELECTROMAGNETIC COMPATIBILITY

This medical device needs special care as far as electromagnetic compatibility is concerned. It must be installed and used in compliance with the information given in the accompanying documents.

Portable and mobile radiocommunication appliances may affect the operation of the thermometer.

MANUFACTURER'S GUIDE AND DECLARATION - ELECTROMAGNETIC EMISSIONS		
The clinical thermometer can be used in the following electromagnetic environment. The customer or user of the clinical thermometer must make sure that the device is used in said environment.		
Emission test	Compliance	Guide to the electromagnetic environment
Irradiated emissions/ CISPR11	Group 1	The clinical thermometer uses RF energy for its internal operation only. Consequently its RF emissions are very low and cause no interference close to electronic appliances.

MANUFACTURER'S GUIDE AND DECLARATION - ELECTROMAGNETIC IMMUNITY		
The clinical thermometer can be used in the following electromagnetic environment. The customer or user of the clinical thermometer must make sure that the device is used in said environment.		
Emission test	Compliance	Guide to the electromagnetic environment
Electrostatic discharges (ESD) IEC/EN61000-4-2	$\pm 6\text{ kV}$ in contact $\pm 8\text{ kV}$ in the air	The floors should be made of wood, cement or baked clay. If the floors are covered in synthetic material, relative humidity should not exceed 30%.
Magnetic field IEC/EN 61000-4-8	3 A/m	The magnetic field should be the typical magnetic field of a commercial or hospital environment.
Irradiated immunity IEC/EN61000-4-3	3 V/m 80 MHz at 2.5 GHz (in appliances that are not life equipment)	

SYMBOLS

: Caution: read the instructions

: application parts type B

NOTE: All the indications, specifications and drawings are based on the latest information available when this manual was printed. They can be subject to variations or changes with no notice required.

This unit has been manufactured in compliance with the European directives 93/42ECC, notably EN 12470-3, EN 60601-1, EN 60601-1-2

DISPOSAL

The device (including its removable parts and accessories) must not be disposed of together with municipal waste at the end of its life, but in compliance with European Directive 2002/96/EC. Since it must be handled separately from household waste, either carry it to a separately-collected waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new device with the same purpose. Any infringement will be severely prosecuted.

The batteries used in this device must be disposed of in the special bins at the end of their life.

Famidoc Technology Co.,Ltd. Add: East 2/F, Zhixiang Building, 71 Block Baoan District, Shenzhen 518101, P.R. China

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Distributed with brand JOYCARE, by JOYCARE S.p.A.,
Via Fabio Massimo, 45 - Roma (Italia)

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GB: WARRANTY CARD

Joycare S.p.A., a company with a sole shareholder, offers a period of 3 year warranty after the date of purchase. The warranty does not cover damage caused by accidental breakage, transportation, incorrect maintenance or cleaning, misuse, tampering or repairs carried out by non authorised third parties. The following is also excluded: damage caused by incorrect product installation, and damage caused by wear and tear, obviously including parts subject to wear and tear. In order to obtain the replacement/repair of products under warranty in cases other than those mentioned above, but which are instead caused by manufacturing defects, please contact your local importer/distributor. Goods that are not held to be in good condition because of the aforesaid causes must be accompanied by this certificate, duly filled-in, as well as a receipt for the goods (invoices and transport documents) and purchasing documentation.

Please keep the original packaging of the product.

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JOYCARE S.p.A. socio unico