BILANCIA BABY

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BABY SCALE

Thanks for having purchased the Joycare Baby Scale, thanks to which you can weigh your baby whenever you want and at home.

Before starting using the unit, it is advisable to read carefully this instruction manual. It is also advisable to save this manual for a future use.

WARNINGS

Use this appliance only for the purposes shown in this instruction manual. Any other purpose should be considered dangerous. Before starting using the appliance, please check that the product is in perfect conditions, without damages. In case of doubts, please address to a specialised centre. Never try to repair the unit by yourself, but always address to a specialised centre. Always place and store the product away from heating sources, sunlight, dampness and sharp objects. When the unit is not used for a long period of time, it is advisable to remove the battery as a possible leaking can damage the product. Always store the unit out of children reach. This appliance may be used by children of more than 8 years of age, persons with physical, sensory or mental capabilities or inexpert persons only if they have been previously trained in using it in a safe manner, and only if they have been informed of the dangers the product can cause. Children may not play with the appliance. Children may not clean or carry out any maintenance on the appliance without adult supervision. This scale is not to be used to determine the weight of objects or substances involved in commercial transactions, for the manufacture of medicines, for calculating tolls, tariffs, taxes, premiums, fines, wages, indemnities or other, similar charges calculated on the basis of weight. This scale can be used to control the weight of children, but it should not be used for diagnostic activities or medical treatments. Weights which appear unusual or abnormal should always be discussed with the child's physician. The baby scale should never be used by adults, but only by children who weight no more than 20 Kg, in order to avoid damaging the product's hardware and/or software, and especially to keep from distorting the calibration, all of which would void the guarantee due to improper use. In order to avoid danger to the user of this device, the doses of any drug or therapy prescribed by one's own doctor should in no circumstances be modified on the basis of data calculated by this device. NB: to avoid possible injury, do not step on the edge of the platform. Packaging components (plastic bags, cardboard, polystyrene etc.) must be kept out of the reach of children as they have the potential to be dangerous and their disposal should be carried out according to the current regulations.

SPECIFICATIONS

- Equipped with a high precision "stain-gauge" sensor and high capability CPU to ensure the weighing accuracy.

- Tare function.
- Auto power off
- Show low voltage
- Capacity: 20KG – Division : 10g
- Unit: KG/LB/ST
- Battery: 2pcs AAA Battery (test batteries are included)
- Work temperature: 0°C- 40°C
- Display usage: LCD Show, 3 units

Operation use conditions:Room temperature:

from 10°C to 40°C. Relative humidity: from 30% to 75% without condensation. Pressure: from 700 – 1060 hPa.

Environmental transportation and storage conditions:

Temperature: from 10 to 40°C. Humidity: from 5 to 95%. Pressure: from 700 – 1060 hPa.

KEY FUNCTION

- ON/OFF/TARE key: Power on, Power off, Power tare
- UNIT key: select the measuring unit for kg/lb/st.

OPERATION

- Take out the scale, and open up the battery cover, install the batteries to the battery case.
- Put the scale on a flat table, and ensure it steady.

WEIGHTING STEP

- 1- Under the state press the key of “ON/OFF/TARE” to power on, the scale will be zero “0.00” appear on the display after 2 seconds.
- 2- Press the UNIT key, you can change “kg”, “lb”, “st”.
- 3- Under the weight state, put the baby on the scale, LCD will show the weight of baby.
- 4- Press the key of “ON/OFF/TARE”, LCD shows “TARE”. Clear the weight of the container, then put the baby into container, the weight will appear on the display.
- 5- Press the key of “ON/OFF/TARE” again , power off.

ERROR DISPLAY INTRODUCTION

When the scale is over loaded, the “EEEE” will appear on the display.

When the screen shows “ LO”, it means that the power isn’t enough, please change the batteries as soon as possible.

ADVICE FOR USE AND CARE

Put the scale on a clean and flat table.

The max. weight of the scale is 20000g, which included the weight before tare.

To ensure the life of the scale, do not put the scale under wet, hot & cold condition.

The product belongs to high precision electronic weighing apparatus, please do not press or drop the scale otherwise the scale is easy to damage.

If the scale surface is dirty, please clean the scale with soft cloth, do not use any chemical detergent and allow water to get inside the scale.

The product is for family weighing only, do not use for legal trade.

If the result appears improperly, please check about whether the batteries and the touching point contact well. Replace for new batteries if used for a long time already, then try again.

MANTENANCE AND CARE

Clean the unit with a wet cloth, avoiding water comes inside the unit. Never use chemical detergents.

Never dip the unit in water. Never place the scale upright as this will cause battery exhaustion. Avoid shaking, dropping or shocking the scale. This is a precision instrument and must be handled with care. Never try to repair the scale by yourself but always address to specialised centres. On the contrary, the warranty will not valid anymore. When the scale is switch off, never place or lay any objects on it.

In case functioning trouble may occur:

- please check if battery has been correctly installed in its compartment
- please check that the unit has been placed on a flat and stable surface.
- please lay your baby on the tray only after the display has shown 0.00

If problems persist, contact the Service Centre JOYCARE.

 This item is made in conformity with all the applicable European directives

DISPOSAL

The device (including its removable parts and accessories) must not be disposed of together with municipal waste at the end of its life, but in compliance with European Directive. Since it must be handled separately from household waste, either carry it to a separately collected waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new device with the same purpose. Any infringement will be severely prosecuted.

The batteries used in this device must be disposed of in the special bins at the end of their life.

Specifications and designs are based on the latest information available at the time of printing and subject to change without notice.

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